

DEBAMED® Absorber with compartments
*Sample tube holder with absorbent properties, for
transporting samples of biological substances.*

DEBATIN®
SAFE. PERFECT. PACKED.



L.E.S.S. FRANCE



► **DEBAMED®** absorber with compartments, *for transporting
samples of biological substances. In order to comply with current ADR and IATA regulations, absorbers must be
positioned between the primary and secondary packaging. The absorber with compartments ensures primary
containers reach their destination safely – even for fragile primary containers (e. g., glass tubes).*

DEBAMED® *Absorber with compartments*
Sample tube holder with absorbent properties, for
transporting samples of biological substances.



Product information

Description

Product characteristics	<ul style="list-style-type: none"> • Ultra-absorbent pad for transporting samples of biological substances in compliance with packaging instructions P650 and P620 • Can be used for all biological fluids (including blood) • Absorbed biological fluids are turned into gel • Protects sample tubes (primary packaging) • No print • Single use only
Standard range	<p>Absorber with 5 compartments</p> <ul style="list-style-type: none"> • Dimensions: 275 x 125 mm • Per compartment: 20 ml blood / 37 ml fluid (e. g.: urine) • Quantity: 400 pcs per box <p>Absorber with 8 compartments</p> <ul style="list-style-type: none"> • Dimensions: 435 x 125 mm • Per compartment: 20 ml blood / 37 ml fluid (e. g.: urine) • Quantity: 250 pcs per box
Country of origin	<ul style="list-style-type: none"> • UK

