

DEBAMED® *Speci-Bag*
Waterproof bag for transporting
samples of biological substances.

DEBATIN
SAFE. PERFECT. PACKED.



L.E.S.S. FRANCE



► **The DEBAMED®** *Speci-Bag* can be used for transporting substances internally around hospitals, or as external secondary packaging for transporting samples of biological substances category B (UN3373) in compliance with the international and European Guidance on regulations for the transport of infectious substances (ADR – WHO – ISO 15189).

Product information

Description

Product characteristics	<ul style="list-style-type: none"> Waterproof transport solution with additional document pouch for forms Easy to handle and quick to open (no risk of injury) Suitable for single-use Tested by an external test lab in accordance with packaging instruction P650 for ADR class 6.2
Material	<ul style="list-style-type: none"> Polyethylene – low environmental impact
Standard range	<ul style="list-style-type: none"> 165 x 265 mm (A5 size), 2,000 pcs per box 245 x 345 mm (A4 size), 1,000 pcs per box 380 x 435 mm (A3 size), 500 pcs per box
Print on standard range	<p>All:</p> <ul style="list-style-type: none"> Biohazard logo <p>Product-specific versions:</p> <ul style="list-style-type: none"> Without print Coloured field Black print Instructions in 4 languages (French, English, German, Spanish)
Closure/opening	<ul style="list-style-type: none"> Strong, waterproof adhesive closure 2 tear notches on the left and right side of the bag (easy to open – no knife or scissors required)
Optional extras	<ul style="list-style-type: none"> Personalised print (logo, name, address, colour) (minimum order: 10,000 pcs) Customised size (minimum order: 10,000 pcs) Absorbent liner (as insert or supplied separately) Without document pouch Supplied in blocks (of 50 or 100 pcs) Bags can be manufactured and printed in a choice of colours Labels with the UN3373 logo printed in accordance with ADR guidelines for outer packaging available.
Country of origin	<ul style="list-style-type: none"> Germany France

